Report Status: Final



Patient Information	Specimen Information	Client Information
DOB: AGE: Gender: Phone: Patient ID:	Specimen: Requisition: Lab Ref #: Collected: Received: Reported:	

COMMENTS: FASTING:NO

Test Name	In Range	Out Of Range	Reference	Range	Lab
POLIOVIRUS (TYPES 1,3) AB					TXC
NEUTRALIZATION					
POLIO 1 TITER	>1:128				
POLIO 3 TITER	>1:128				
REFERENCE RANGE:					
POLIO 1 TITER:	<1:8				
POLIO 3 TITER:	<1:8				

The presence of neutralizing serum antibodies (titers of 1:8 up to >1:128) against polioviruses implies lifelong immunity. Some persons without detectable titers (<1:8) may also be immune as demonstrated by elicitation of a secondary-type serum antibody response upon rechallenge with live polio vaccine.

This test was developed and its analytical performance characteristics have been determined by Focus Diagnostics. It has not been cleared or approved by FDA. This assay has been validated pursuant to the CLIA regulations and is used for clinical purposes.

PERFORMING SITE:

TXC FOCUS DIAGNOSTICS, 33608 ORTEGA HIGHWAY BLD B-WEST WING, SAN JUAN CAPISTRANO, CA 92675-2042 Laboratory Director: HOLLIS BATTERMAN, MD, CLIA: 05D0644251