

Patient Information	Specimen Informa	ntion	Client Information	
Test Name	In Range	Out Of Range	Reference Range	Lab
MYASTHENIA GRAVIS PANEL 3 ACETYLCHOLINE RECEPTOR				EZ
BINDING ANTIBODY	<0.30		nmol/L	112
Reference Ranges for Acetyl	choline Recep	tor		
Binding Antibody:				
Negative: < or =0.30 nmol/	L			
Equivocal: 0.31-0.49 nmol/ Positive: > or =0.50 nmol/				
ACETYLCHOLINE RECEPTOR				ΕZ
BLOCKING ANTIBODY ACETYLCHOLINE RECEPTOR	<15		<15 % Inhibition	ΕZ
MODULATING ANTIBODY	11	_	% Inhibition	
		ference Range: 2% INHIBITION	•	
mled a secret of 3 3 3 3		1		
This test was developed and characteristics have been d			CS.	
It has not been cleared or been validated pursuant to				
for clinical purposes.		Tacions and is t	iseu	
STRIATED MUSCLE AB SCREEN	NEGATIVE	V	NEGATIVE	EZ
This test was developed and				
characteristics have been d It has not been cleared or				
been validated pursuant to				
for clinical purposes. STRIATED MUSCLE AB TITER				EZ

Test Not Performed. Screening test Negative or Not Detected. Titer not performed.