

Patient Information	Specimen Information	n Client Information
DOB: AGE: Gender: Phone: Patient ID:	Specimen: Requisition: Lab Ref #: Collected: Received: Reported:	
Test Name CALPROTECTIN, STOOL		Out Of RangeReference RangeLab683.9 Hmcg/g

<15.625 - 50 mcg/g Normal >50 - 120 mcg/g Borderline >120 mcg/g Abnormal

Calprotectin in Crohn's disease and ulcerative colitis can be five to several thousand times above the reference population (50 mcg/g or less). Levels are usually 50 mcg/g or less in healthy patients and with irritable bowel syndrome. Repeat testing in 4-6 weeks is suggested for borderline values.

## **PERFORMING SITE:**

Quest, Quest Diagnostics, the associated logo and all associated Quest Diagnostics marks are the trademarks of Quest Diagnostics.