

Specimen ID:
Control ID:

Acct #: 17452095 Phone: (800) 539-6119 Rte: 00
Walk-In Lab, LLC
1645 Tiffany Lane
Mandeville LA 70448

Patient Details

DOB:
Age(y/m/d):
Gender: SSN:
Patient ID:

Specimen Details

Date collected:
Date entered:
Date reported:

Physician Details

Ordering: D ELWOOD
Referring:
ID: 1033316112
NPI: 1033316112

General Comments & Additional Information

Alternate Control Number:
Total Volume: Not Provided

Alternate Patient ID: Not Provided
Fasting: No

Ordered Items

proBNP; Blood Drawing

TESTS	RESULT	FLAG	UNITS	REFERENCE	INTERVAL	LAB
proBNP						
proBNP	788	High	pg/mL	0 - 210		01
Please Note:						01

The following cut-points have been suggested for the use of proBNP for the diagnostic evaluation of heart failure (HF) in patients with acute dyspnea:

Modality	Age (years)	Optimal Cut Point
Diagnosis (rule in HF)	<50	450 pg/mL
	50 - 75	900 pg/mL
	>75	1800 pg/mL
Exclusion (rule out HF)	Age independent	300 pg/mL

01	BN	LabCorp Burlington 1447 York Court, Burlington, NC 27215-3361	Dir: William F Hancock, MD
02	MB	LabCorp Birmingham 1801 First Avenue South, Birmingham, AL 35233-1935	Dir: Brian Ragland, MD

For inquiries, the physician may contact **Branch: 504-828-2666 Lab: 800-762-4344**